











SURGI STEEL

Multifilament Stainless Steel




Wire Suture B.P.

NEEDLE			LENGTH	GAUGE
TEMPORARY PACEMAKER WIRES - STANDARD				32 swg
16mm				
Curved round bodied			45cm	1904D
63mm				
Straight cutting			if break point required	B1904D
25mm				
1/2 circle round bodied			45cm	1920
25mm				
1/2 circle round bodied			45cm	1924D
90mm				
Curved cutting			if break point required	B1924D

32 swg soft braided, specially insulated stainless steel wire. 7.5cm uninsulated at either end. 30cm of insulation. 45cm overall length. The 63mm needle used can be supplied with a 'breakpoint' if required. 'Breakpoint' needle pacemaker wires are supplied on 59.5cm wire (52 cm of insulation and 7.5cm uninsulated at the non-'breakpoint' end).

SURGI STEEL

Monofilament Stainless Steel Wire Suture B.P.

NEEDLE	LENGTH	GAUGE	
MONOFILAMENT STEEL WIRE SUTURE		1	5
		(4)	(7)
40mm			
1/2 circle cutting point - heavy Multipack (4 sutures per envelope) 		45cm	1714M
55mm			
1/2 circle Trocar point - heavy 		75cm	1798H