

MONOSORB - Monofilament Polydioxanone Suture

w Monofilament Polydioxanone E.P. Violet.

w Monosorb suture is a synthetic absorbable sterile surgical suture prepared from the polyester poly (dioxanone).

w Monosorb sutures are dyed with D&C Violet # 2.

w Monosorb is available in a range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.

w Minimal tissue reaction.

w Easy passage through tissue.

w Uses - Polydioxanone sutures are intended for use where an absorbable suture or liouature is indicated.

w Sterilised by Ethylene Oxide.

w Monosorb sutures from Sutures Ltd. are guaranteed to meet the general requirements of the European Pharmacopoeia (Ph.Eur.).

SURGICAL INDICATIONS suggested uses	NEEDLE	LENGTH	GAUGE																	
			8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (4.5)	1 (6)	2 (8)								
	8mm																			
VASCULAR	1/4 circle spatula double	45cm				DW10379D														
	Curved round bodied double	45cm				DW10669D														
	12mm																			
	Curved cutting - fine	45cm				DW10219	DW10220													
	13mm																			
VASC/PAED/UROL	Curved round bodied	45cm				DW10386														
	Curved round bodied	70cm					DW10387	DW10388S												
	15mm																			
	Curved cutting - fine	45cm					DW10250	DW10251												
	16mm																			
UROL/VASC	1/2 circle round bodied - double	90cm					DW10800D	DW10801D												
UROL/G.I./PAED	Curved round bodied	70cm					DW10790	DW10791	DW10792											
SKIN/ SUBCUTICULAR	Curved Cutting	70cm					DW10860	DW10861												
	19mm																			
SKIN/ SUBCUTICULAR	Curved reverse cutting	45cm						DW10901R	DW10902R											
	20mm																			
PAED/G.I./UROL	1/2 circle round bodied	70cm						DW10921	DW10922	DW10923										
	Curved cutting - fine	45cm					DW10265	DW10266												
	22mm																			
ORAL/DENTAL	1/2 circle cutting	45cm								DW10577										
	25mm																			
PLASTIC/SKIN/ SUBCUTICULAR	Curved cutting - fine	45cm						DW10971												
	26mm																			
	1/2 circle round bodied	70cm						DW10941	DW10942	DW10943	DW10944									
SKIN/ SUBCUTICULAR	Curved reverse cutting	45cm							DW10832R											
	30mm																			
G.I./PAED/FASCIA	1/2 circle round bodied	70cm						DW10951	DW10952	DW10953	DW10954	DW10955								
	35mm																			
	1/2 circle round bodied	70cm							DW10473	DW10474	DW10475									
	Curved reverse cutting	70cm							DW10563R	DW10564R										
	40mm																			
G.I./PAED/FASCIA	1/2 circle round bodied	70cm							DW10512	DW10513	DW10514	DW10515S7								
G.I./PAED/FASCIA	1/2 circle round bodied	90cm									DW10514S8	DW10515								
FASCIA	1/2 circle reverse cutting	90cm								DW10523R	DW10524R	DW10525R								
	45mm																			
	1/2 circle round bodied	90cm																	DW10855	
	50mm																			
FASCIA/ PERITONEUM	1/2 circle round bodied	90cm							DW10063	DW10064	DW10065	DW10066								
FASCIA	1/2 circle reverse cutting	90cm								DW10544R	DW10545R									
	60mm																			
SUBCUTICULAR	Straight cutting	70cm						DW10571	DW10572	DW10573										
NON-NEEDED MONOSORB																				
	Cut lengths	1.5m								DW10502	DW10503	DW10504	DW10505							

KEY: VASC = VASCULAR, PAED = PAEDIATRIC, UROL = UROLOGICAL

KEY: VASC = VASCULAR, UROL = UROLOGICAL

KEY: UROL = UROLOGICAL, G.I. = GASTRO INTESTINAL, PAED = PAEDIATRICAL

KEY: UROL = UROLOGICAL, G.I. = GASTRO INTESTINAL, PAED = PAEDIATRICAL

KEY: G.I. = GASTRO INTESTINAL, PAED = PAEDIATRICAL

KEY: G.I. = GASTRO INTESTINAL, PAED = PAEDIATRICAL

KEY: G.I. = GASTRO INTESTINAL, PAED = PAEDIATRICAL