

PLAIN CATGUT - PLAIN CATGUT B.P.

- w Catgut Plain B.P.
- w Manufactured from uniform, twisted strands of collagen.
- w Packaged in a hydrating fluid containing isooctanol and water.
- w Origin - Purified animal intestines (bovine).
- w Type - Absorbable by enzymatic digestion.
- w Uses - Plain catgut sutures are intended for use where an absorbable suture or ligature is indicated.
- w Plain catgut is available in a range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.
- w Sterilised by Gamma irradiation.
- w Plain catgut sutures from Sutures Ltd. are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur.) and the United States Pharmacopoeia (U.S.P.).

SURGICAL INDICATIONS suggested uses		NEEDLE	LENGTH	GAUGE												
				7/0	6/0	5/0	4/0	3/0	2/0	0	1	2	3	4	5	6
				(0.2)	(0.17)	(0.13)	(0.2)	(0.25)	(0.33)	(0.4)	(0.5)	(0.6)	(0.7)	(0.8)	(0.9)	(1.0)
8mm																
STRABISMUS	1/4 circle spatula - double		45cm			DW159D	DW160D	DW161D								
10mm																
PAEDIATRIC/PLASTIC/ ORTHOPAEDIC	Curved cutting - double		40cm				DW180D	DW181D								
CONJUNCTIVAL	Curved reverse cutting		45cm	DW188R	DW189R											
STRABISMUS	Curved reverse cutting - double		45cm			DW219RD										
13mm																
PAEDIATRIC/SKIN/PLASTIC	Curved cutting		45cm				DW230	DW231								
16mm																
PAEDIATRIC/UROLOGICAL	Curved round bodied		75cm				DW190+ DW250	DW251	DW252	DW253	DW254					
PAEDIATRIC/UROLOGICAL	Curved round bodied - double		75cm					DW251D								
PAEDIATRIC/UROLOGICAL	1/2 circle round bodied		75cm					DW261	DW262	DW263						
PAEDIATRIC/PLASTIC/ SUBCUTICULAR/SKIN	Curved cutting		75cm					DW801	DW802	DW803	DW804					
STRABISMUS	Curved cutting - double		38cm					DW281D	DW282D							
DENTAL/SKIN	1/2 circle cutting		45cm					DW236	DW237							
DENTAL/SKIN	1/2 circle cutting		75cm					DW291	DW292							
19mm																
SKIN/PAEDIATRIC	Curved reverse cutting		50cm				DW790R	DW791R	DW792R							
SKIN/PAEDIATRIC	Curved cutting		45cm						DW792	DW793						
20mm																
G.I./UROLOGICAL	1/2 circle round bodied		75cm					DW361	DW362	DW363						
22mm																
SKIN	Curved cutting		75cm						DW372							
SKIN/DENTAL	1/2 circle cutting		75cm							DW383	DW384					
25mm																
UROLOGICAL/PAEDIATRIC	1/2 circle round bodied		75cm				DW390	DW391	DW392	DW393	DW394	DW395				
DENTAL	1/2 circle cutting point		75cm							DW833	DW834					
28mm																
SKIN/DENTAL	Curved cutting		75cm				DW820	DW821	DW822	DW823	DW824	DW825	DW826			
SKIN/DENTAL	Curved reverse cutting		45cm						DW822R	DW823R						
30mm																
UROLOGICAL/GENERAL USE	1/2 circle round bodied		75cm					DW441	DW442	DW443	DW444	DW445	DW446			
GENERAL CLOSURE/ ORTHOPAEDIC/ VETERINARY	1/2 circle round bodied - heavy		75cm											DW445H		
SUBCUTICULAR/ PAEDIATRIC	Curved cutting		75cm						DW462	DW463	DW464					
35mm																
UROLOGICAL/ GENERAL CLOSURE	1/2 circle round bodied		75cm						DW232	DW233	DW234	DW235				
GENERAL CLOSURE/ HENNA/OBS & GYN/ ORTHOPAEDIC/ VETERINARY	1/2 circle round bodied - heavy		90cm							DW239H	DW234H	DW235H	DW236H			
GENERAL CLOSURE/ OBS & GYN/UROLOGICAL	1/2 circle cutting point		75cm						DW492	DW493	DW494	DW495	DW496			
40mm																
SUBCUTICULAR/ GENERAL CLOSURE/ OBS & GYN	1/2 circle Trocar point		75cm							DW613	DW941					
SUBCUTICULAR/ GENERAL CLOSURE/ OBS & GYN	1/2 circle round bodied		75cm						DW622	DW623	DW624	DW625	DW626			
GENERAL CLOSURE/ OBS & GYN/ORTHOPAEDIC	1/2 circle round bodied - heavy		90cm							DW623H	DW624H	DW625H	DW626H			
SUBCUTICULAR/OBS & GYN/ VETERINARY	1/2 circle reverse cutting		75cm									DW645R	DW646R			
50mm																
G.I./GENERAL CLOSURE/ FAT STITCH	1/2 circle round bodied		75cm							DW663						
GENERAL CLOSURE/ FAT STITCH	1/2 circle round bodied - heavy		75cm							DW663H	DW664	DW665	DW666			
60mm																
SKIN	Straight cutting		75cm								DW683					
NON-NEEDED PLAIN CATGUT																
LIGATION/ THREADING EYED NEEDLES	Ligature lengths		5 x 38cm								DW33	DW34				
LIGATION/ THREADING EYED NEEDLES	Cut lengths		75cm				DW60	DW61	DW62	DW63	DW64	DW65	DW66			
LIGATION/ THREADING EYED NEEDLES	Standard lengths		1.5m				DW100	DW101	DW102	DW103	DW104	DW105	DW106			