

SURGIBOND - Coated Braided Polyester Suture

w Central Resiloid Inductor Gamma R.P. (Diam)

w Tose - Non-absorbable synthetic material.

w Surgibond is available in a range of tissue sizes and lengths, non-needled or attached to stainless steel needles of various types and sizes.

w Manufactured from a smooth, coated polyester material with a constant curve in the Sutures 1 (needle)

w Extremely strong for lino term support.

w Soft and stable providing good knot security.

w Minimal tissue reaction.

w Easy passage through tissue and reduced tissue trauma.

w Uses - Cardiac, vascular and orthopaedic procedures where lino term support from a braided non-absorbable suture is essential.

w Sterilized by Gamma Irradiation.

w Surgibond Sutures from Sutures Ltd. are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur.) and the United States Pharmacopoeia (U.S.P.).

SPECIAL INSTRUCTIONS	NEEDLE	LENGTH	GAUGE														
			0/0	0/0	7/0	6/0	5/0	4/0	3/0	2/0	1/0	0	0	0			
			18	16	14	12	10	9	8	7	6	5	4	3	2	1	0
VASCULAR	Carved round bodied - double	45cm	DW71030	DW71060		DW71100											
PLASTIC/ORTHOPAEDIC	1/2 circle reverse cutting	45cm		DW71288													
TENDON REPAIR	1/2 circle cutouts - double	45cm			DW71290	DW71350	DW71360										
RETINAL	1/4 circle spatula - double	45cm				DW71300	DW71310										
VASCULAR	Carved round bodied	75cm				DW7149											
DAGVICYSTORHINOSTOMY	5/8 circle cutting point - double	white 25cm					DW71500W										
PLASTIC	Carved reverse cutting	45cm					DW71598										
PLASTIC/SUBCUTICULAR/TENDONS	Carved cutting	45cm				DW7709	DW7710	DW7711									
VASCULAR/PAEDIATRIC/TENDONS	Carved round bodied	75cm				DW7039	DW7040	DW7041									
VASCULAR/PAEDIATRIC/TENDONS	Carved round bodied - double	75cm				DW70300	DW71900	DW70410									
VASCULAR/PL/PAEDIATRIC/TENDONS	1/2 circle round bodied	75cm					DW7060	DW7061	DW7062	DW7063							
VASCULAR/PL/PAEDIATRIC/TENDONS	1/2 circle round bodied - double	60cm					DW72300	DW72310	DW72320								
VASCULAR/PL/PAEDIATRIC/TENDONS	1/2 circle round bodied - double	75cm				DW70500	DW70600	DW70610	DW70620								
CARDIO/VASCULAR	1/2 circle round bodied - double	90cm							DW72300	DW72340							
VASCULAR	Carved round bodied	75cm					DW7070	DW7071									
VASCULAR	Carved round bodied - double	75cm					DW70200	DW70210	DW70220								
PLASTIC/SUBCUTICULAR	Carved cutting	45cm					DW7220	DW7221	DW7222								
VASCULAR	1/2 circle cutting point	90cm				DW7229	DW7300	DW7301	DW7302	DW7303							
CARDIO/VASCULAR	1/2 circle cutting point - double	90cm				DW72900	DW73000	DW73010	DW73020	DW73030							
PLASTIC	Carved cutting	45cm					DW7230	DW7231	DW7232	DW7233							
VASCULAR	1/2 circle round bodied	75cm					DW7311	DW7312	DW7313	DW7314							
CARDIO/VASCULAR/PL	1/2 circle round bodied - double	75cm					DW73300	DW73310	DW73320	DW73330							
VASCULAR	1/2 circle cutting point	75cm					DW7341	DW7342	DW7343								
CARDIO/VASCULAR	1/2 circle cutting point - double	75cm					DW7410	DW7420	DW7430								
ORTHOPAEDIC	1/2 circle trapez point	75cm					DW7061	DW7062									
VASCULAR	1/2 circle trapez point - double	75cm					DW76610	DW76620									
PL/CAERIAL/VASCULAR	1/2 circle round bodied	75cm					DW7381	DW7382	DW7383	DW7384							
PL/CAERIAL/VASCULAR	1/2 circle round bodied - heavy	1m					DW73810	DW73820	DW73830	DW73840							
ORTHOPAEDIC	1/2 circle round bodied - heavy	75cm					DW7391	DW7392	DW7393	DW7394							
PLASTIC	Carved cutting	45cm					DW7420	DW7421	DW7422								
PLASTIC	Carved cutting	75cm					DW7391	DW7392	DW7393	DW7394							
PLASTIC	Carved cutting - fine	45cm					DW7391	DW7392	DW7393	DW7394							
TENDONS/VASCULAR	Carved cutting point	25cm					DW7401	DW7402	DW7403	DW7404							
VASCULAR	1/2 circle cutting point	75cm					DW7411	DW7412	DW7413	DW7414							
VASCULAR	1/2 circle cutting point - double	75cm					DW74110	DW74120	DW74130	DW74140							
VASCULAR	1/2 circle cutting point - fine	45cm					DW7431	DW7432	DW7433								
VASCULAR	1/2 circle cutting point - fine - dual	90cm					DW74310										
PLASTIC	Carved cutting	45cm					DW7271	DW7272	DW7273								
PLASTIC	Carved reverse cutting	75cm					DW7420										
GENERAL CLOSURE/HENNA/ GYN REPAIR/ORTHOPAEDIC	1/2 circle round bodied	75cm					DW7443	DW7443	DW7444	DW7445							
HENNA/GYN REPAIR/ORTHOPAEDIC	1/2 circle round bodied - double	1.3m					DW7450										
VASCULAR/CAERIAL/ORTHOPAEDIC	1/2 circle trapez point	75cm					DW7672	DW7673	DW7674								DW7676
VASCULAR/CAERIAL/ORTHOPAEDIC/GYN	1/2 circle cutting point	75cm					DW7703	DW7703	DW7704	DW7705	DW7706						
VASCULAR/CAERIAL	1/2 circle cutting point - double	75cm								DW77050	DW77050						
VASCULAR/CAERIAL	1/2 circle cutting point - double	90cm					DW7751	DW7752	DW7753	DW7754	DW7755						
PLASTIC	Carved cutting	75cm					DW7751	DW7752	DW7753	DW7754	DW7755						
GENERAL CLOSURE/HENNA/ GYN REPAIR/ORTHOPAEDIC	1/2 circle round bodied	75cm					DW7472	DW7473	DW7474	DW7475	DW7476						
VASCULAR	1/2 circle round bodied - double	75cm					DW7470	DW7471	DW7472	DW7473	DW7474						
HENNA/GYN REPAIR/ORTHOPAEDIC	1/2 circle round bodied - heavy	75cm					DW7470H	DW7471H	DW7472H	DW7473H	DW7474H						
VASCULAR/CAERIAL	1/2 circle cutting point - double	90cm					DW7480	DW7481	DW7482	DW7483	DW7484						
ORTHOPAEDIC/GYN REPAIR	1/2 circle cutting point - heavy	75cm					DW7480H	DW7481H	DW7482H	DW7483H	DW7484H						
SKIN CLOSURE/ DRAIN & LINE FIXATION/ VASCULAR	Carved cutting point	75cm					DW7603	DW7604	DW7605	DW7606							
GENERAL CLOSURE/HENNA	1/2 circle round bodied	1m					DW7503			DW7505							
ORTHOPAEDIC/GYN REPAIR/ HENNA/GENERAL CLOSURE/ VETERINARY	1/2 circle round bodied - heavy	75cm								DW7504H	DW7505H	DW7506H					
ORTHOPAEDIC/ GENERAL CLOSURE/ VETERINARY	1/2 circle cutting point - heavy	75cm									DW7534H	DW7535H	DW7536H				
GENERAL CLOSURE/ VETERINARY	1/2 circle trapez point	75cm									DW7684	DW7685	DW7686				
GENERAL CLOSURE	1/2 circle cutting point	75cm									DW7554						
SKIN CLOSURE	Straight cutting	90cm									DW7573						
SKIN CLOSURE	Straight cutting	75cm									DW7593	DW7593					
SKIN CLOSURE	Straight cutting	1m									DW7521						
WOUND RETENTION	Carved cutting - double	75cm															DW7644
WOUND RETENTION	Carved cutting + vinyl tubing Multipack - 2 per packet 12 packets per box	50cm															DW7645
SKIN CLOSURE/ DRAIN & LINE FIXATION	Carved cutting	1m									DW7643	DW7644					

NON-NEEDED SURGIBOND

USGATION/ THREADING EYED NEEDLES	Cut lengths	13 X 60cm								DW7011	DW7012	DW7013	DW7014				
USGATION/ THREADING EYED NEEDLES	Cut lengths	1.8m								DW7001	DW7002	DW7003	DW7004	DW7005	DW7006		
USGATION/ THREADING EYED NEEDLES	Cut lengths	2m								DW7021	DW7022	DW7023	DW7024	DW7025	DW7026		