

SURGILACTIN - Coated Braided Polyglactin 910 B.P. Violet and Undyed

w Coated Braided Polyglactin 910 B.P. Violet and Undyed.

w Surgilactin suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide.

w Surgilactin sutures are dyed with D&C Violet # 2.

w Surgilactin is available in a range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.

w Minimal tissue reaction.

w Easy passage through tissue.

w Uses - Surgilactin sutures are intended for use where an absorbable suture or fixation is indicated.

w Sterilized by Ethylene Oxide.

w Coating - Surgilactin sutures are coated with a mixture of copolymer of glycolide & lactide and Calcium Sterate.

w Surgilactin sutures from Sutures Ltd. are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph.Eur.), and the United States Pharmacopoeia (U.S.P.) *except for an occasional slight overage in some names.*

SURGICAL INDICATIONS expected use	NEEDLE	LENGTH	GAUGE											
			8/0 (9.4)	7/0 (10.5)	6/0 (12.7)	5/0 (15)	4/0 (17.5)	3/0 (20)	2/0 (23)	0 (25.6)	1 (28)	2 (31)	3 (34)	
8mm														
GLAUCOMA/CATARACT	Curved spatula - double	30cm	DW116270	DW116280										
	Curved ultra spatula - double	30cm	DW116377	DW116380										
8mm														
GLAUCOMA/CONJUNCTIVAL	Curved round bodied	30cm	DW116675	DW11668										
GLAUCOMA/CATARACT	1/2 circle spatula - double	45cm		DW113480	DW113490									
STRABISMUS	1/4 circle spatula - double	45cm			DW113790	DW113800								
	1.2mm													
	Curved reverse cutting	45cm				DW11220R								
	Curved reverse cutting	75cm					DW11881R							
	1.8mm													
UBOL /PAED/GI	Curved round bodied	60cm		DW113865G	DW113875G	DW11388								
	1.6mm													
UBOL /PAED/GI	Curved round bodied	45cm					DW11791							
SKIN/ SUBCUTICULAR	Curved cutting	75cm				DW11860	DW11861	DW11862	DW11862					
UBOL /PAED/GI	1/2 circle round bodied	75cm				DW11800	DW11801	DW11802	DW11802					
	Curved reverse cutting	75cm				DW11860R	DW11861R	DW11862R						
	1.8mm													
SKIN/ SUBCUTICULAR	Curved reverse cutting	45cm						DW11901RS						
	Curved reverse cutting	70cm							DW11902R	DW11903R				
	Curved reverse cutting	75cm						DW11901R						
	20mm													
UBOL	1/2 circle round bodied	75cm					DW11921	DW11922	DW11923					
	24mm													
	Curved reverse cutting	70cm					DW11936R	DW11937R	DW11938R					
	26mm													
SKIN/ SUBCUTICULAR	Curved reverse cutting	45cm							DW11832RS					
	Curved reverse cutting	50cm					DW11831R							
	1/2 circle reverse cutting	70cm						DW11832R						
	Curved cutting	75cm						DW11832	DW11833	DW11834				
SKIN/ SUBCUTICULAR	1/2 circle round bodied	75cm					DW11941	DW11942	DW11943					
	30mm													
OBS/GYN/FACIA	1/2 circle cut taper	75cm							DT11703	DT11704	DT11705			
	Curved reverse cutting	70cm							DW11813R					
	1/2 circle round bodied	70cm												
G.I./PAED	1/2 circle round bodied - heavy	75cm					DW11951	DW11952	DW11953	DW11954	DW11955			
	1/2 circle round bodied - heavy	75cm								DW11960	DW11961			
	35mm													
	Curved cutting	75cm							DW11562	DW11563				
OBS/GYN/FACIA	1/2 circle cut taper	75cm							DT11873					
OBS/GYN/FACIA	1/2 circle cut taper - heavy	75cm							DT11873H		DT11565	DT11566		
G.I./PAED	1/2 circle round bodied	75cm							DW11993	DW11994	DW11995			
	38mm													
	1/2 circle round bodied	75cm							DW11882					
	40mm													
	1/2 circle reverse cutting	70cm											DW11526RS	
	1/2 circle cutting	70cm											DW11526S	
FACIA/ORTHO	1/2 circle cut taper	75cm							DT11415	DT11416	DT11417	DT11418		
G.I./FACIA	1/2 circle round bodied	75cm							DW11512	DW11513	DW11514	DW11515	DW11516	
FACIA/ PERTONEUM	1/2 circle reverse cutting	75cm								DW11524R	DW11525R			
	1/2 circle mitting	75cm								DW11523	DW11524	DW11525		
OBS/GYN/G.I.	1/2 circle trocar point heavy	75cm										DW11825		
	45mm													
FACIA/ORTHO	1/2 circle cut taper heavy	75cm								DT11533H	DT11536H			
SKIN/ SUBCUTICULAR	Curved reverse cutting	75cm							DW11583R					
OBS/GYN/FACIA	1/2 circle round bodied heavy	75cm								DW11854	DW11855	DW11856		
	50mm													
FACIA/ PERTONEUM	1/2 circle round bodied	75cm							DW11553	DW11554	DW11555	DW11556		
	60mm													
SKIN/ SUBCUTICULAR	Straight cutting	75cm							DW11572	DW11573				

NON-NEEDED SURGILACTIN

Cut lengths	17 x 45cm						DW11401	DW11402						
Cut lengths	10 x 50cm						DW11461	DW11462	DW11463	DW11464				
Cut lengths	1.5m						DW11501	DW11502	DW11503	DW11504	DW11505	DW11506		

KEY: UROL = UROLOGICAL, PAED = PAEDIATRIC, G.I. = GASTRO INTESTINAL

KEY: UROL = UROLOGICAL, PAED = PAEDIATRIC, G.I. = GASTRO INTESTINAL

KEY: UROL = UROLOGICAL, PAED = PAEDIATRIC, G.I. = GASTRO INTESTINAL

KEY: UROL = UROLOGICAL

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL

KEY: G.I. = GASTRO INTESTINAL, PAED = PAEDIATRIC

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL

KEY: G.I. = GASTRO INTESTINAL, PAED = PAEDIATRIC

KEY: G.I. = GASTRO INTESTINAL

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL, GASTRO INTESTINAL

KEY: ORTHO = ORTHOPAEDIC

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL