

SURGIQUICK - Coated Braided Polyglycolic acid B.P. Undyed and Violet

w Braided Polyglycolic acid B.P. Violet and Undyed.

w Surgiquick suture is a synthetic rapidly absorbable sterile surgical suture composed of homopolymer of glycolide.

w Surgiquick sutures are dyed with D&C Violet # 2.

w Surgiquick is available in a range of gauge sizes and lengths, non-neededled or attached to stainless steel needles of varying types and sizes.

w Minimal tissue reaction.

w Easy passage through tissue.

w Uses - Surgiquick sutures are intended for use where an absorbable suture is indicated for short term support.

w Sterilised by Ethylene Oxide.

w Coating - Polyglycolic acid is a homopolymer of glycolic acid, which is braided into Sutures of various sizes. The braid is coated with a mixture of Calcium Stearate and Polycaprolactone.

w Surgiquick sutures from Sutures Ltd. are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph.Eur.), and the United States Pharmacopoeia (U.S.P.).

SURGICAL INDICATIONS suggested uses	NEEDLE	LENGTH	GAUGE																		
			8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)									
	12mm																				
SKIN/SUBCUTICULAR	P Needle curved reverse cutting	45cm					DP12220														
	13mm																				
	Straight cutting	45cm							DW12398												
	16mm																				
	Bi-curved reverse cutting	45cm						DW12004R	DW12005R												
	Curved cutting	75cm						DW12860	DW12861	DW12862											
	19mm																				
SKIN/SUBCUTICULAR	P Needle curved reverse cutting	75cm							DP12901R	DP12902R											
SKIN/SUBCUTICULAR	Curved reverse cutting	75cm							DW12901R	DW12902R											
	22mm																				
ORAL/DENTAL	1/2 circle cutting	45cm							DW12576	DW12577											
	26mm																				
SKIN/SUBCUTICULAR	P Needle curved reverse cutting	75cm								DP12832RS7	DP12833R										
SKIN/SUBCUTICULAR	Curved reverse cutting	75cm								DW12832RSU7	DW12833RU										
	35mm																				
SKIN/SUBCUTICULAR	P Needle curved reverse cutting	75cm								DP12562R	DP12563R										
SKIN/SUBCUTICULAR	Curved reverse cutting	75cm								DW12562R	DW12563RU										
EPISIOTOMY	1/2 circle cut taper	90cm									DT12873US9	DT12874S9	DT12565S								
	40mm																				
	1/2 circle cut taper	90cm										DT12415S	DT12416S	DT12417S							
	45mm																				
	1/2 circle cut taper	90cm										DT12533S9	DT12534S								
	Curved reverse cutting	75cm										DW12583R									
	Straight cutting	45cm							DW12007												
	60mm																				
	Straight cutting	75cm								DW12572	DW12573										
NON-NEEDED SURGIQUICK																					
	Cut lengths	1.5m							DW12500	DW12501	DW12502	DW12503	DW12504	DW12505							