

SURGISORB - Coated Braided polyglycolic acid B.P. violet

- w Braided Polyglycolic acid B.P. Violet and Undyed.
- w Surgisorb suture is a synthetic absorbable sterile surgical suture composed of homopolymer of glycolide.
- w Surgisorb sutures are dyed with D&C Violet # 2.
- w Surgisorb is available in a range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.
- w Minimal tissue reaction.
- w Easy passage through tissue.
- w Uses - Surgisorb sutures are intended for use where an absorbable suture or ligature is indicated.
- w Sterilised by Ethylene Oxide.
- w Coating - Polyglycolic acid is a homopolymer of glycolic acid, which is braided into Sutures of various sizes and coated for smooth passage through tissue and easy knot tie down.
- w Surgisorb sutures from Sutures Ltd. are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph.Eur.), and the United States Pharmacopoeia (U.S.P.).

SURGICAL INDICATIONS <small>synongeted uses</small>	NEEDLE	LENGTH	GAUGE															
			8/0 (0.9)	7/0 (0.8)	6/0 (0.7)	5/0 (1)	4/0 (1.3)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)						
GLAUCOMA/CATARACT	Curved spatula - double	30cm			DW8628D													
	Curved ultra spatula - double	30cm	DW8637D	DW8638D														
	8mm																	
	Curved round bodied	30cm	DW8667S															
	Curved round bodied	45cm	DW8667	DW8668	DW8669													
	1/2 circle spatula - double	45cm	DW8348D		DW8349D	DW8350D												
	1/4 circle spatula - double	45cm			DW8379D	DW8380D	DW8381D											
	12mm																	
STRABISMUS	Curved reverse cutting	45cm				DW8220R	DW8881RS											
	Curved reverse cutting	75cm					DW8881R											
	13mm																	
IRCOL/PAED/G.I.	Curved round bodied	60cm					DW8388											
IRCOL/PAED/G.I.	Curved round bodied	45cm				DW8790	DW8791											
SKIN/SUBCUTICULAR	Curved cutting	75cm		DW8859	DW8860	DW8861	DW8862	DW8863										
IRCOL/PAED/G.I.	1/2 circle round bodied	75cm			DW8800	DW8801	DW8802	DW8803										
	Curved reverse cutting	75cm				DW8860R	DW8861R	DW8862R	DW8863R									
SKIN/SUBCUTICULAR	13mm																	
	Curved reverse cutting	45cm					DW8901RS											
	Curved reverse cutting	70cm						DW8902R	DW8903R									
	Curved reverse cutting	75cm					DW8901R	DW8902RS7										
	20mm																	
IRCOL	1/2 circle round bodied	75cm					DW8921	DW8922	DW8923									
24mm																		
	Curved reverse cutting	70cm				DW8935R	DW8936R	DW8937R	DW8938R	DW8939R								
PLASTIC/SUBCUTICULAR	26mm																	
	Curved reverse cutting	45cm							DW8832RS									
	Curved reverse cutting	50cm					DW8831R											
	1/2 circle reverse cutting	70cm						DW8852R	DW8853R									
SKIN/SUBCUTICULAR	Curved cutting	75cm						DW8832	DW8833	DW8834								
	1/2 circle round bodied	75cm					DW8941	DW8942	DW8943		DW8944	DW8945						
	30mm																	
OBS/GYN/FASCIA	1/2 circle cutting point	75cm						DW8702	DW8703	DW8704	DW8705							
	Curved reverse cutting	70cm						DW8813R										
	1/2 circle round bodied	70cm								DW8954S						DW8956S7		
G.I./PAED	1/2 circle round bodied	75cm				DW8951	DW8952	DW8953	DW8954	DW8955	DW8956							
OBS/GYN/FASCIA	1/2 circle round bodied heavy	75cm							DW8960	DW8961								
PLASTIC/SKIN/SUBCUTICULAR	35mm																	
	Curved cutting	75cm						DW8562	DW8563									
	1/2 circle cutting point	75cm							DW8873	DW8874	DW8875							
OBS/GYN/FASCIA	1/2 circle cutting point heavy	75cm							DW8564	DW8565	DW8566							
OBS/GYN/FASCIA	1/2 circle round bodied	75cm						DW8992	DW8993	DW8994	DW8995							
G.I./PAED	38mm																	
	1/2 circle round bodied	75cm						DW8882	DW8883									
	40mm																	
ORTHOGYN	1/2 circle reverse cutting	70cm															DW8526BS	
	1/2 circle cutting	70cm															DW8526S	
	1/2 circle cutting point heavy	75cm							DW8415	DW8416	DW8417							
G.I./FASCIA	1/2 circle round bodied	75cm						DW8513	DW8514	DW8515	DW8516							
FASCIA/PERITONEUM	1/2 circle reverse cutting	75cm							DW8548	DW8549	DW8550						DW8526R	
OBS/GYN/G.I.	1/2 circle cutting	75cm						DW8522	DW8523	DW8524	DW8525	DW8526						
	1/2 circle trocar point heavy	75cm									DW8825							
	45mm																	
FASCIA/ORTHO	1/2 circle cutting point heavy	75cm						DW8533H	DW8534H	DW8535H	DW8536H							
SKIN/SUBCUTICULAR	Curved reverse cutting	75cm						DW8583R										
OBS/FASCIA	1/2 circle round bodied heavy	75cm							DW8854	DW8855	DW8856							
FASCIA/PERITONEUM	50mm																	
	1/2 circle round bodied	75cm							DW8553	DW8554	DW8555	DW8556						
	60mm																	
SKIN/SUBCUTICULAR	Straight cutting	75cm					DW8570		DW8572	DW8573								
NON-NEELED SURGISORB																		
	Cut lengths	17 x 45cm						DW8401	DW8402	DW8403								
	Cut lengths	10 x 50cm						DW8461	DW8462	DW8463	DW8464							
	Cut lengths	1.5m						DW8500	DW8501	DW8502	DW8503	DW8504	DW8505	DW8506				