

SURGICAPRONE - Monofilament Polyglecaprone 25 B.P. Violet and Undyed

w Monofilament Polyglecaprone 25 B.P. Violet and Undyed.

w Suricaprone suture is a synthetic absorbable sterile surgical suture composed of Poly (polycolide-co-ε-caprolactone).

w Suricaprone sutures are dyed with D&C Violet # 2.

w Suricaprone is available in a range of gauge sizes and lengths, non-needled or attached to stainless steel needles of varying types and sizes.

w Minimal tissue reaction.

w Easy passage through tissue.

w Uses - Suricaprone sutures are intended for use where an absorbable suture or ligature is indicated for short term support.

w Sterilised by Ethylene Oxide.

w Suricaprone sutures from Sutures Ltd. are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.) and the European Pharmacopoeia (Ph.Eur.).

SURGICAL INDICATIONS suggested uses	NEEDLE	LENGTH	GAUGE									
			6/0 (6.7)	5/0 (1)	4/0 (4.5)	3/0 (3)	2/0 (1.5)	1/0 (1)	1/2 (0.7)			
		13mm										
PLASTIC/ SUBCUTICULAR	P Needle curved reverse cutting	45cm			DP13002							
	Curved reverse cutting	45cm	DW13759RU	DW13760RU								
	1/2 circle round bodied	45cm	DW13753									
		16mm										
PLASTIC/ SUBCUTICULAR	P Needle curved reverse cutting	45cm			DP13007	DP13008						
	Curved reverse cutting	45cm			DW13860RSU	DW13861RSU						
	1/2 circle round bodied	70cm										
PAED/UROL/G.I. ENT/UROL	Curved round bodied	70cm				DW13013						
	1/2 circle round bodied	70cm				DW13018	DW13019					
	19mm											
PLASTIC/ SUBCUTICULAR	P Needle curved reverse cutting	45cm				DP13023	DP13024					
	Curved reverse cutting	45cm			DW13900RSU	DW13901RU	DW13902RU					
	1/2 circle round bodied	70cm				DW13028	DW13029	DW13030				
	25mm											
PAED/MUSCLE	1/2 circle cut taper	70cm									DT13033	
	26mm											
PLASTIC/ SUBCUTICULAR	P Needle curved reverse cutting	70cm					DP13037					
	Curved reverse cutting	70cm				DW13831RSU7	DW13832RSU7	DW13833RS				
	1/2 circle round bodied	70cm					DW13942	DW13943				
	1/2 circle round bodied - double	70cm					DW13042D	DW13043D				
	27mm											
	5/8 circle round bodied	70cm					DW13047					
G.I.	1/2 circle round bodied - blunt violet	70cm								DW13953BS		
	32mm											
G.I./PAED	1/2 circle round bodied	70cm					DW13057	DW13058	DW13059	DW13060		
	35mm											
MUSCLE/PAED/G.I.	Curved round bodied	70cm						DW13063	DW13064			
	35mm											
G.I.	1/2 circle round bodied - blunt	70cm						DW13068	DW13069			
	1/2 circle round bodied - blunt heavy	90cm									DW13075	
	36mm											
OBS/GYN	1/2 circle round bodied - blunt heavy	90cm										DW13075
	1/2 circle cut taper	70cm						DT13078	DT13079	DT13080		
GYN/MUSCLE	36mm											
	1/2 circle round bodied	70cm						DW13083	DW13084	DW13085		
	1/2 circle reverse cutting	70cm						DW13088R	DW13089R			
G.I./MUSCLE	40mm											
	1/2 circle round bodied - blunt	90cm							DW13089	DW13100		
	1/2 circle round bodied - blunt heavy	90cm									DW13103	
	1/2 circle reverse cutting - heavy	90cm									DW13109R	
MUSCLE	1/2 circle round bodied - heavy	70cm								DW13114	DW13115	
	45mm											
MUSCLE/OBS/GYN	1/2 circle round bodied - blunt heavy	90cm										DW13120
	1/2 circle round bodied - heavy	90cm										DW13125
MUSCLE/OBS/GYN	50mm											
	1/2 circle round bodied - blunt heavy	90cm										DW13130
MUSCLE/GYN	1/2 circle round bodied - heavy	90cm										
	60mm											
SUBCUTICULAR	Straight cutting	70cm				DW13571	DW13572	DW13573				
	Straight cutting with 2 beads & 2 collars	70cm					DW13572X	DW13573X				
	70mm											
	1/2 circle round bodied - blunt	90cm										DW13140

KEY: PAED = PAEDIATRIC, UROL = UROLOGICAL, G.I. = GASTRO INTESTINAL
KEY: UROL = UROLOGICAL

KEY: PAED = PAEDIATRIC

KEY: G.I. = GASTRO INTESTINAL

KEY: G.I. = GASTRO INTESTINAL, PAED = PAEDIATRIC

KEY: G.I. = GASTRO INTESTINAL, PAED = PAEDIATRIC

KEY: G.I. = GASTRO INTESTINAL

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL

KEY: GYN = GYNAECOLOGICAL

KEY: G.I. = GASTRO INTESTINAL

KEY: GYN = GYNAECOLOGICAL, ORTHO = ORTHOPAEDIC

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL

KEY: OBS = OBSTETRICAL, GYN = GYNAECOLOGICAL

KEY: GYN = GYNAECOLOGICAL

NON-NEEDED SURGICAPRONE

cut lengths	1.5m				DW13502	DW13503	DW13504	DW13505
-------------	------	--	--	--	---------	---------	---------	---------